		Twist						Twist R			
ration	Possible set-ups	Rear mounting	Rear mounting with Caregiver Controller	Front mounting	Front mounting with Handlebar Control	Front mounting with Caregiver Controller	Rear mounting	Rear mounting with Caregiver Controller			
Configu	Accessories needed for the set-up	Frame: code KT-16-0000 + Controller right side: code KT-CTR-DX-0000 or Controller left side: code KT-CTR-SX-0000	Frame: code KT-16-0000 + Caregiver Controller: code KT-PA-0000	Controller right side: code KT-CTR-DX-0000 or Controller left side: code KT-CTR-SX-0000	Handlebar Control: code KT-25-0000	Caregiver Controller: code KT-PA-32-0000	Controller right side: code KT-CTR-DX-0000 or Controller left side: code KT-CTR-SX-0000	Caregiver Controller: code KT-PA-32-0000			
	Compatible with rigid and folding wheelchairs with 24-inch wheels	YES (*2)	YES (*2)	YES (*2)	YES (*2)	YES (*2)	YE	S (*2)			
iit	Compatible with rigid and folding wheelchairs with 26-inch wheels	NO	NO	YES (*2)	YES (*2)	YES (*2)	YES (*2) NO				
atib	Compatible with wheelchairs with wheels smaller than 24-inch wheel	NO	NO	YES (*2)	YES (*2)	YES (*2)					
튭	Compatible with tilt-in-space wheelchairs	Only available for 24-inch wheels (*2)	Only available for 24-inch wheels (*2)	YES / available for all wheels (*2)	NO	YES (*2)	Possible with 24	or 26-inch wheels			
ن	Compatible with wheelchairs with an inner wheel-to-wheel distance	From 330 to 510 mm (*3)	From 330 to 510 mm (*3)	All sizes	All sizes	All sizes	From 330 to	510 mm (*3)			
	Weight with battery + accessories installed	9,4 kg	9,4 kg	8,7 kg	11,2 kg	8,7 kg	7,6 kg				
	Weight without battery - device only			7,3 kg (only Twist)			6,2 kg				
	Dimensions			290 x 355 x 155 mm			290 x 350 x 155 mm				
	Max. travel range (*1)		Max 20 km (*1)								
	Approved for a maximum rated slope		Maxin		Maximum rated slope 10% and maximum user weight of 120 Kg						
	Maximum negotiable slope (*4)	15%	15%	10%	15%	15%	15%	15%			
ies	Maximum speed	10 km/h	6 km/h	6 km/h	20 km/h	6 km/h	10 km/h	6 km/h			
alit	Maximum speed in reverse	1	1	1	5 km/h	1	1	1			
i i	Electronic brake adjustable from Web App Klaxon®	NO	NO	YES	YES	NO	1	10			
Ē	Electronic brake at the accelerator release adjustable or deactivable via Web App Klaxon®	NO	NO	NO	YES	YES	NO	YES			
and	Mechanical disc brake	NO	NO	NO	YES	NO	NO				
ıres	Mechanical parking brake	NO	NO	NO	YES	NO	!	10			
eatı	Start & Stop function via Controller	YES	NO	YES	NO	NO	YES	NO			
_	Emergency electronic brake for immediate motor stop	NO	YES	YES	NO	YES	NO	NO			
	Emergency propulsion shutdown via Controller	YES	NO	YES	NO	NO	YES	NO			
	Speed selection in 4 levels via Controller	YES	NO	YES	NO	NO	YES	NO			
	Speed control via accelerator lever	NO	YES	NO	YES	YES	NO	YES			
	Downhill speed limiter and energy recover	NO	NO	YES	YES	NO	NO	NO			
	Automatic maintenance of set speed when going uphill	YES	NO	YES	NO	NO	YES	NO			
NS NS	Automatic speed reduction when cornering	NO	NO	YES	NO	NO	NO	NO			
afet	Automatic emergency shutdown in the event of tipping	YES	YES	YES	YES	YES	YES	YES			
s s	Mechanical directionality control system	NO	NO	YES	NO	NO	NO	NO			
_	Battery	Lithium-ion 36V /144Wh / 4A									
nica Ita	Motor				700W						
ech D <sub>a</sub>	Wheel	10 inch with replaceable tyre									
_	Connettivity	Bluetooth / WiFi									
tion	Twist's product certification	. \									
tifica	Klaxon Mobility GmbH Company's certification										
Cert	Software's Certification	"IEC 62304 (Medical device software - Software life cycle processes)"									

		Twist						Twist R	
		Rear mounting	Rear mounting with Caregiver Controller	Front mounting	Front mounting with Handlebar Control	Front mounting with Caregiver Controller	Rear mounting	Rear mounting with Caregiver Controller	
	Indoor driving	••••	••••	••••	•••••	••••	•••••	••••	
	Driving in urban environments with a smooth surface	••••	••••	••••	••••	••••	••••	••••	
f us	Driving in urban environments with uneven pavements / cobblestones / slabs	•••	••••	••••	••••	••••	•••	•••	
ns o	Driving on uphill slopes not exceeding 10%	••••	••••	•••	••••	••••	•••••	••••	
tatio	Driving on uphill slopes not exceeding 15%	••••	••••	Not allowed	••••	•••	••••	•••••	
Ē	Driving in presence of light counter slopes		•••••	••••	•••••	••••	•••	••••	
and	Driving in presence of accentuated counter-slopes	•	••••	••••	•••••	••••	•	••••	
ons	Driving on downhill slopes of no more than 10%		•••••	••••	•••••	••••	•••	••••	
dati	Driving on downhill slopes of more than 10%	•	•••	••••	•••••	••••	•	•••	
E E	Driving on off-road tracks such as grass or beaten and even ground	••	•••••	••••	•••••	•••••	••	••••	
E CO	Driving in wet and rainy conditions with no deep puddles present	••	•••••	••		••••	•••	••••	
2	Driving on off-road tracks with gravel, mud, water puddles, stones	Not allowed	Not allowed	Not allowed	Not allowed	Not allowed	Not allowed	Not allowed	
	Driving in wet and rainy conditions with deep puddles present	Not allowed	Not allowed	Not allowed	Not allowed	Not allowed	Not allowed	Not allowed	